

Benchtop Tube Sealer WeldOn-B



REMI Bench top sealer **WeldOn-B** is designed to meet high standards in blood safety. **WeldOn-B** helps to keep the blood bag system completely sterile. It ensures quality sealing with clear rupture line for easy tear-apart segments. Consistent performance, safety features & monitoring system makes it ideally suited for easy & long lasting usage

WeldOn-B

Salient Features

- Senses the blood bag tube for automatic sealing process
- Anti-spark and overheat protection system to avoid damage to system & tubing
- Visual indication for : -
 - Ready to seal,
 - Seal process
 - Power supply,
- Easy to remove sealing head cap for quick cleaning
- Seal with clear break line for snap-apart separation
- Splash guard for technician safety



Seal with clear break line for snap-apart separation

Technical Specifications

| Parameter | WeldOn-B |
|------------------------------|--|
| Dimension(W x D x H) (mm) | 220 X 200 X 120 |
| Weight (Kg) | 4 |
| RF Frequency | 40.68 MHZ |
| RF Output | 80 watt |
| Seal time | 1.5 Sec (auto adjustable) |
| Sealing Type | Automatic |
| EEC Standard | Class 1 (93/42 and 47/2007 MDD) |
| Input Voltage | 115/230 V, 50/60 Hz |
| Ref Standard | E 60601-1 for electric safety, E 60601-1-2 for EMC |

Specifications subject to change without notice

