



# VDRL Shaker VDS 300



**VDS 300**

VDS 300 is specially designed for smooth & uniform mixing of samples, Helps to standardize the process with utmost user convenience

### Salient Features:

- Universal platform with adjustable roller to fix the sample in micro-plate, petri-dish, glass vessels or test tubes
- Brush less motor for maintenance free operations
- Adjustable rollers to fix samples of different sizes
- Rubber cushion help in stable sample holding

### Applications :

- VDRL testing
- Microbiology
- Molecular biology
- Pathology
- Biotechnology

### Technical Specifications

Parameter	VDS 300
Dimension (W x D x H)	270 x 260 x 125 (mm)
Platform size	7" x 11"
Shaking orbit	10 mm
Countdown timer	0-99 minutes
Min speed	50 - 300 RPM
Power Supply	220-240 v, 50 Hz, Single phase