



## REMI Water Bath

### RWB Series



RWB 6

#### Technical Specifications

Parameter	RWB 6
Temperature range	RT $\pm$ 5 °C ~99°C
Temperature stability	$\pm$ 0.5°C
Tracking alarm	$\pm$ 2°C
Chamber volume	15 Ltrs
Interior Dimension (W×D×H, mm)	450×300×110
Timing range	1~999min
Water Bath Cover	Double lines six holes
Electrical requirements	220V, 50Hz, Single phase
Power Consumption	1500W

REMI Water Bath is multipurpose water bath with adjustable double line six holes to hold different size samples for easy & hygienic heating. RWB provides precise and constant temperature helps achieving uniform heating for the samples.

#### Salient Features:

- Microprocessor control with timing function
- Digital display
- Audible and visual alarm for over temperature