



## REMI Microplate Shaker

The Digital Microplate Shaker is an essential instrument for mixing individual micro plates of different type

### Salient Features:

- Tiny footprint-as small as the plate itself
- Digital Display with Speed & Time setting
- Speed selection from 300 to 1800 RPM  
Maintenance free BLDC motor
- Orbital Diameter of 2 mm for thorough mixing, microprocessor controlled for high accuracy
- Silent & stable operation
- Timer setting of to 99 min & Continuous mode
- A Special turn down toggle clamp for spill proof & simple side in/side out for loading / unloading of microplate
- Universal fit-comes with interchangeable clamps for different micro plates (Tool free and easy clamp change procedure)
- Even mixing in all wells across the plate
- Comes With attachment for 20 x 1.5/2ml microtubes
- Can be used inside fume hoods / cold rooms



MPS 1

### Technical Specifications

Specifications	Microplate Shaker
Speed setting	300 to 1800 RPM
Orbital diameter	2 mm
Motor type	Brushless DC motor
Timer setting	Up to 99 mins + infinite mode
Display	Digital display (LED) of Speed & Time
Power supply	Universal (100 to 230V, 50/60 Hz)
Dimensions (W x D x H)	127 x 85 x 86 mm
Weight	1.4 Kg
Standard attachments	<ol style="list-style-type: none"> <li>1. 20 x 1.5/2.0 ml Microtubes (foam)</li> <li>2. Replaceable clamps to use different types of standard micro plates</li> </ol>