



REMI Water Bath RWB Series



RWB 6

Technical Specifications

REMI Water Bath is multipurpose water bath with adjustable double line six holes to hold different size samples for easy & hygienic heating. RWB provides precise and constant temperature helps achieving uniform heating for the samples.

Salient Features:

- Microprocessor control with timing function
- Digital display
- Audible and visual alarm for over temperature

Parameter	RWB 6
Temperature range	RT \pm 5 $^{\circ}$ C ~99 $^{\circ}$ C
Temperature stability	\pm 0.5 $^{\circ}$ C
Tracking alarm	\pm 2 $^{\circ}$ C
Chamber volume	15 Ltrs
Interior Dimension (WxDxH, mm)	450x300x110
Timing range	1~999min
Water Bath Cover	Double lines six holes
Electrical requirements	220V, 50Hz, Single phase
Power Consumption	1500W