



## Portable Donor Chair



**PDC 1X**



Folded

REMI Portable Donor Chair model PDC 1X, specially designed for blood donation camps to provide comfort and safety to the blood donors.

- Ergonomic robust design
- Stable & aesthetically beautiful
- Complete construct with stainless steel
- High quality anti bacterial fabric
- Foldable for easy transport
- Adjustable safety belt for patient safety
- Towel bar at top
- 15 kg light weight for easy transportation
- IV Stand optional (standalone)
- Comfortable seating
- Trendelenburg (Head low) & Reverse-Trendelenburg (Sitting Posture) position
- Elegant ergonomic knobs for the hand rest adjustment
- Adjustable head rest positions for easy phlebotomy
- Elegant ergonomic knob for hand rest adjustment
- Adjustable hand rest positions (Up / down and horizontal / circular)
- Comfortable molded foamed hand rests

### Technical Specification

Model	PDC 1X
Seat width	600 mm
Length	1700 mm
Seat Height from floor	1300 mm
Weight Capacity	150 kg
Hand rest dimension (W x L mm)	150 x 450
Net Weight	15 kg