



Plasma Expressor **XPRESS - M**



XPRESS - M is specially designed plasma expressor for separating out the blood components. The front panel of the

XPRESS - M is loaded with high tension spring for quick transfer of plasma into other bag.

XPRESS - M is ergonomically designed for long lasting trouble free operations

Salient Features

- User friendly with see through front panel to allow visibility of layer movement
- High tension spring loaded thick compression plate for fast plasma expression
- Ergonomic handle with auto adjustable lever for easy operation

Technical Specifications

Parameter	Unit	Xpress M
Dimensions (W x D x H)	mm	165x230x280
Weight	kg	3.5
Body Construct		Powder Coated CRCA Steel
Front Panel		Acrylic