



## Automated Plasma Expresser Xpress A



**Xpress A**

Xpress A is a semi-automated plasma expressor for quick & easy component separation. The optical sensing of RBC's & automatic clamping reduces manual error & avoids RBC contamination.

### Salient Feature

- Automated fast acting clamp mechanism, no RBC contamination
- Spring loaded transparent compression plate for quick plasma expression
- Adjustable spring tension
- Ergonomic handle for easy operation
- LED power indication, buzzer alarm on clamping
- Suction holders for firm base grip
- Flexibility of manual operation
- Stainless steel framing available (Optional)

### Technical Specifications

Overall Dimensions	167x270x290mm (LXBXH)
Max Tube Diameter	6.5mm
Weight	4 kg
Operating Temperature	0 to 40°C
Input Supply	110 – 240v