



Therapeutic Chairs



Comfort

Hygiene

Safety



Therapeutic Chairs



1. Head rest adjustable for any comfortable position
2. Cushioned body supports for extra comfort
3. Movable (horizontal & vertical) fully padded broad armrests with glass holders
4. Adjustable / motorized foot rest (Dzire 4X) with removable transparent plastic cover
5. Anti-fungal, Anti-bacterial, easy to clean leatherette upholstery
6. Sturdy & stable steel construction, allows load carrying capacity up to 200 kgs.
7. Best in class Dual lock, Anti-skid, Anti-static PU castors
8. Tissue paper / Towel holder
9. 2/3/4 independent motors for height, back, leg & Foot movements
10. Infusion stand, adjustable both sides
11. Emergency switch for head low position



Therapeutic Chairs

The innovative design of the Dzire 2X / 3X / 4X Therapeutic Chair underlines the high quality standards. The comfortable upholstery creates a pleasant environment and feeling of wellbeing for patients. Two/ Three/ Four electric motors offer simple and smooth movements of all desired positions further adding to patient relaxation whilst receiving treatment.

The remote control provides greater flexibility as all desired chair positions can be adjusted easily.

Dzire Xtreme Comfort

- Specially designed for comfortable physiological posture
- Non Deformable Soft PU foam
- Dual padded Back Cushion
- Adjustable Head Support
- Vertical & Horizontal Armrest Movement
- Ergonomically designed for Comfortable Posture
- Height adjustment –user friendly for clinical team (Dzire 3X / 4X)
- Smooth back rest and foot rest movement
- Easily movable infusion stand


Dzire Xtreme Hygiene

- Removable cushions for cleaning & maintenance
- Epoxy powder coated robust frame with rounded corners
- Anti-fungal, Anti-bacterial, easy to clean upholstery
- Elegant, easy to clean ABS moulded cover to prevent dust accumulation on mechanism



Dzire Xtreme Safety

- Single touch button for shock position - ensures patient safety
- Anti-static, Anti-skid, Dual lock PU castors
- Sturdy, stable & non-tilt design capable to take weight up to 200 kg





Therapeutic Chairs



Technical Specifications

Parameters	Dzire 2X	Dzire 3X	Dzire 4X
No of Motors / Actuators	2	3	4
Motorized movement	Leg rest, Back rest	Leg rest, Back rest, Height	Leg rest, Back rest, Height, Foot rest
Foot rest	Spring loaded & self adjustable with removable cover		Motorized with removable cover
Wheels	Anti-static, Anti-skid, Dual lock PU castors		
Construction	Epoxy powder coated steel frame		
Removable cushion	Anti-fungal, Anti-bacterial, leatherette with CFC free PU foam		
Attachments	Infusion stand, adjustable both sides		
Hand rest movement	Vertical & Horizontal movement and Height adjustable		
Weight carrying capacity	200 kg		
Weight	75kg	80kg	85kg
Compliance	IEC 60601-1, IEC 60601-1-2		
Electrical	220v-240v 50Hz, Single phase		

* Battery Backup (Optional)

