



MS 800 Plus

Ultra - flat compact magnetic stirrer, guaranteed with modern magnet coil technology. Wear - free with no moving parts. To ensure better mixing, the lab disc can reverse direction of rotation automatically every 30 seconds

- No moving parts - truly maintenance free
- Low aspect ratio (just 12 mm thick)
- Software driven gradual acceleration
- Adjustable speed (15 to 1500 RPM) with digital display
- Microprocessor controlled with last run memory feature
- Stirring cap. 800 ml - set up place of 100 mm
- Pulse mode - reverses the direction every 30 secs for proper mixing
- Universal power supply - small wall mount adaptor
- IP 65 compliant design (Water spill resistant)
- Corrosion resistant SS top surface to fight harsh chemicals
- Recommended Stir bar length 25 mm
- Supplied with 1 (one) high strength stir bar

| Technical Data | |
|--------------------------------------|---|
| Specifications | MS-800 Plus |
| Max RPM | 1500 RPM |
| Speed Setting | 15 to 1500 RPM (at every 50 RPM) |
| Set - up plate material | Stainless Steel |
| Speed indicator | Yes |
| Permissible ambient temperature [°C] | 5 - 40 |
| Permissible relative humidity [%] | 80 |
| Load capacity (Water) | 800 ml |
| Recommended Stir bar length | 25 mm |
| Pulse function | Yes |
| Sync Mode | No |
| Last run memory | Yes |
| Dimensions (W X D X H) | 120 x 164 x 12 mm |
| Weight | 400 gm |
| Protection Class | IP 65 compliant design (Water spill resistant) |

Supply: 220-240 Volts, 50 Hz Single Phase.