

Model	Stirring Shaft		Impeller		Motor	Speed (rpm)		Stirring capacity		Digital Speed Indicator
	Length (mm)	Ø (mm)	Type	Ø (mm)		Min	Max	Type of fluid	Max. Volume	
Direct Drive Stirrers										
RQ-121/D	250	8	Axial Turbine	36	PMDC	200	2000	Water/Light Oil	10 Ltrs.	✓
RQ-126/D	450	10	Pitched Fan	75	PMDC	120	2000	Water/Light Oil	25 Ltrs.	✓
RQ-129/D	550	10	Pitched Fan	100	PMDC	100	1400	Water/Light Oil	50 Ltrs.	✓
RQM-122/R	250	6	Pitched Fan	38	AC/DC 1/20 HP	800	4000	Water	7 Ltrs.	✗
RQ-122/D	250	6	Pitched Fan	38	AC/DC 1/20 HP	800	4000	Water	7 Ltrs.	✗
RQ-124A/D	350	8	Pitched Fan	52	AC/DC 1/8 HP	300	4000	Water	20 Ltrs.	✗
RQT-124A/D	350	8	Pitched Fan	52	AC/DC 1/8 HP	300	4000	Water	20 Ltrs.	✓
RQ-134S/D	450	12	Pitched Fan x 2	55	AC/DC 1/4 HP	900	5000	Water	35 Ltrs.	✗
RQ-134L/D	450	12	Pitched Fan x 2	65	AC/DC 1/2 HP	900	6000	Water	60 Ltrs.	✗
RQ-134H/D	450	12	Pitched Fan x 2	75	AC/DC 2/3 HP	500	5000	Water	80 Ltrs.	✗
RQT-134H/D	450	12	Pitched Fan x 2	75	AC/DC 2/3 HP	500	5000	Water	80 Ltrs.	✓
Geared Stirrers										
RQG-121/D	250	6	Axial Turbine	36	PMDC	50	550	Light Oil	5 Ltrs.	✓
RQG-126/D	450	8	4 Cross Blade	75	PMDC	50	750	Paint/Gear Oil	15 Ltrs.	✓
RQG-129/D	550	10	4 Cross Blade	100	PMDC	50	550	Paint/Gear Oil	25 Ltrs.	✓
RQ-128/D	310	8	4 cross Blade	90	AC/DC 1/20 HP	50	150	Light Oil	3 Kgs.	✗
RQ-128B/D	310	8	4 cross Blade	90	AC/DC 1/8 HP	50	150	Light Oil	7 Kgs.	✗
RQ-128D/D	450	12	4 Cross Blade x 2	120	AC/DC 1/2 HP	50	280	Light Oil	20 Kgs.	✗
RQG-128A/D	450	8	4 Cross Blade	75	AC Brushless	50	550	Water	10 Ltrs.	✗
RQG-128C/D	450	10	4 Cross Blade	90	AC Brushless	50	500	Paint/Gear Oil	20 Ltrs.	✗
AC Stirrers / Flameproof Stirrers										
RQ-140/D	730	10	4 Cross Blade	100	AC Brushless	100	1400	Water/Light Oil	40 Ltrs.	✓
RQ-140/DFP	550	10	4 Cross Blade	100	Flame Proof	100	1400	Water/Light Oil	40 Ltrs.	✓
Homogenisers / Emulsifiers										
RQ-127/D	250	8	Saw Cutter	30	AC/DC 1/8 HP	300	6000	Water	8 Ltrs.	✗
RQT-127/D	250	8	Saw Cutter	30	AC/DC 1/8 HP	300	6000	Water	8 Ltrs.	✓
RQ-127A/D	275	8	Teflon Pestle		AC/DC 1/8 HP	300	8000		upto 100ml	✗
RQT-127A/D	275	8	Teflon Pestle		AC/DC 1/8 HP	300	8000		upto 100ml	✓
RQ-130/D	480	12	Saw Cutter	40	AC/DC 1/4 HP	300	5000	Water	20 Ltrs.	✗
RQT-130/D	480	12	Saw Cutter	40	AC/DC 1/4 HP	300	5000	Water	20 Ltrs.	✓
RQ-130B/D	480	12	Saw Cutter	60	AC/DC 2/3 HP	300	5000	Water	35 Ltrs.	✗
RQT-130B/D	480	12	Saw Cutter	60	AC/DC 2/3 HP	300	5000	Water	35 Ltrs.	✓
RQ-140/DE	480	12	Saw Cutter	40	AC Brushless	500	5000	Water	35 Ltrs.	✓

Supply: 220-240 Volts, 50 Hz Single Phase.